


Magnetic Nano Particles for In-Vitro Diagnostics

Ferrotec Corporation has utilized its many years of experience in creating a “**Magnetic Nano-particles Developer Kit**” for in-vitro biomedical applications. The kit offers a unique combination of surface modified and unmodified superparamagnetic particles to researchers and engineers for study of microbead synthesis, hyperthermia, control drug delivery and MRI contrast agent. In the kit are included four types of surfactant coated dry particles with dispersability in various solvents, two types of aqueous based ferrofluids and one bare particles suspension.

This application note describes the main features of the Developer kit, solvents for different products and provides a bibliography of published work.



TEM picture of EMG1200

Dry surfactant coated particles when dissolved in an appropriate solvent form a stable



Ferrofluid is attracted to a magnet as one homogeneous phase


colloid. These and other products in the kit (except for EMG1111) are referred to as ferrofluids. A short description of ferrofluids follows.

A ferrofluid is a stable suspension of nano order size magnetic particles in various carrier liquids. These particles are coated with stabilizers to prevent agglomeration. The physical properties of ferrofluids, such as volatility, environmental capability and viscosity, etc., are mostly determined by the choice of the carrier.

When a magnetic field is applied, a ferrofluid acquires a net magnetic moment due to the orientation of particles in the field direction. These fluids exhibit no remanence or hysteresis. The maximum obtainable moment per unit volume is proportional to the loading density of the magnetic material and is characterized in ferrofluid products as the saturation magnetization value. At saturation, all the particles have their magnetic moments aligned with the external field. Above saturation, the magnetization of fluid is independent of the applied field.


The magnetic particles contained in ferrofluids have a nominal diameter of 10 nm (1/100000 of 1 mm) and are single domain. Some distribution of particles may, however, be noted in TEM measurements. Stable colloidal suspensions in water are available from 1.7 volume percent to 5.0 volume percent of magnetic particles. Also higher concentrated ferrofluids for some particles in hydrocarbon solvent can be prepared. The cores of these particles are made of iron oxides which are known to be compatible with living tissues. From the X-ray measurements, the iron oxides are identified to be a mixture of Fe_3O_4 and $\gamma\text{-Fe}_2\text{O}_3$.

In biomedical applications, the magnetic particles in ferrofluids are used either directly or as a component of polymeric supports for biologically active species. The magnetic properties of these supports allow them to be rapidly separated from complex liquid mixtures even by relatively weak magnetic fields. Depending on the nature of the biologically active species absorbed on the surface of the support; viruses, bacteria, and various cells may be selectively separated.




Magnetization curve for typical iron oxide nano particles. no remanence or hysteresis

Therapeutic agents have been incorporated into or onto the magnetically active polymer particles. The particles are concentrated by magnetic fields at specific body sites where they deliver high local concentration of drug. Ferrofluids have also been used to incorporate magnetic particles into preformed biologically active polymer gels which are then used as supports for affinity chromatography. The simplified separation is made possible due to the magnetic properties of the gels and eliminates the usual centrifugation and column chromatography steps.



Particle size distribution for typical iron oxide, counted from TEM picture

The advantage of using magnetic polymers in affinity chromatography lies in that they allow a rapid absorption of molecules from colloidal solutions or those containing cell debris. Conventional purification methods do not allow for such a quick retrieval. This method is particularly valuable if the requirement is to rapidly isolate labile enzymes or enzyme complexes after their liberation from the cell. It is likely that once the enzyme is affinity-absorbed within the pores of a gel, it is more protected, for example, against proteolysis enzymes.



XRD data for typical iron oxide particles

Ferrotec Corporation offers a starter **“Magnetic Nano-particles Developer Kit”**

consisting of (1) a selection of dry particles having different types of surface modifications, (2) water base ferrofluids stabilized with anionic or cationic surface active agents and (3) a slurry of magnetic particles. The slurry serves two main purposes: first, it allows the researcher to coat the particles with his/her own surfactant and second, some researchers require uncoated particles as reference sample. Nominal properties of each product in the kit are summarized in Tables I.



Appearance of “Magnetic Nano Particles Developer Kit”

Ferrofluid Type	EMG 1200	EMG 1300	EMG 1400	EMG 1500	EMG 607	EMG 707	EMG 1111*
Form	Dry particles				Water base ferrofluid		Water base slurry *
Quantity	2 g	2 g	2 g	2 g	12 ml	12 ml	12 ml
coating nature	fatty acid	polymeric	hydrophobic	polar	Cationic	Anionic	no coating
Typical Particle Concentration in ferrofluid by weight (incl. surfactant)	-	-	-	-	19 %	17 %	19 %
Content of iron oxide in dry particles	about 70 %	about 70 %	about 80 %	about 75 %	-	-	-
Saturation Magnetization of ferrofluid	-	-	-	-	10 mT	10 mT	(about 16 mT)
Typical Saturation Magnetization of particles	82 $\mu\text{Am}^2/\text{kg}$ (60 emu/g)				57 $\mu\text{Am}^2 \text{kg}^{-1}$ (45 emu/g)	57 $\mu\text{Am}^2 \text{kg}^{-1}$ (45 emu/g)	82 $\mu\text{Am}^2 \text{kg}^{-1}$ (65 emu/g)
Initial Susceptibility	0.20	0.22	0.23	0.25	0.36	0.36	0.65
Average Particle	10 nm						
pH	-	-	-	-	9-10	8-9	5-7
Density	-	-	-	-	1100 kg m^{-3}	1100 kg m^{-3}	1180 kg m^{-3}

Table1: Nominal Properties of nano magnetic particles and ferrofluids in the Developer Kit, *EMG 1111 contains aggregates of particles.

The process of dissolving the dry particles in a compatible solvent is simple. Take a small amount of dry particles in a dish or vial and add a solvent of choice from Table II. The particles can be easily dispersed in the solvent by stirring or shaking if the affinity between the particles and solvent is good. If the particles are hard to disperse, heating and stirring by appropriate means may be effective.

EMG 1111 is bare iron oxide particles with no surface treatment, so we recommend an ultra sound treatment (sonication) prior to use of the product.

	water	Methanol	IPA	Butanol	MEK	Methyl acetate	Toluene	Heptane	Xylene
EMG1200	-	-	-	-	-	-	OK	OK	OK
EMG1300	-	-	-	-	-	-	OK	OK	OK
EMG1400	-	-	-	-	-	-	OK	OK	OK
EMG1500	-	-	-	OK	-	-	-	-	-
EMG606	OK	-	-	-	-	-	-	-	-
EMG707	OK	-	-	-	-	-	-	-	-
EMG1111	-	-	-	-	-	-	-	-	-

Table II: Solubility of nano magnetic particles in different solvents

BIBLIOGRAPHY

1. Anderson, L. and Mosbach, K.: Magnetic Ferrofluids for Preparation of Magnetic Polymers and their Application in Affinity Chromatography, Nature, 270, 259 (1977).
2. Czerlinski, G. Senyei, A., and Widder, K.; Magnetic Guidance of Drug Carrying Microspheres, Journal of Applied Physics, 49, (6), 3578, (1978).
3. Giaever, I.; Magnetic Separation of Biological Particles, U.S. Patent Number 3,970,518, July 20, 1976.
4. Hunter, JA.; Ferrography-A New Method for Isolation of Particles from Biological Fluids, Journal of Clinical Pathology, 35,689 (1982).
5. Margel, S., Rembaum, A., and Zisblatt, S.;Polyglutaraldehyde: A New Reagent For Coupling Proteins to Microspheres and For Labeling Cell Surface Receptors. II. Simplified Labeling Method By Means of Nonmagnetic and Magnetic Polyglutaraldehyde Microspheres, Journal of Immunological Methods, 28, 341 (1979).
6. Molday, R.S., Rembaum, A., and Yen, S.P.S.; Application of Magnetic Microspheres in Labeling and Separation of Cells, Nature, 268, 437 (1977).
7. Morimoto, Y.; Magnetic Guidance of Ferro-colloid Entrapped Emulsion for Site Specific Drug Delivery. The Chemical and Pharmaceutical Bulletin, (Tokyo) 1,279 (1983).

8. Morimoto, Y.; Biomedical Applications of Magnetic Fluids II., Journal of Phannacobiodynamics, 4 (8), 624 (1981).
9. Mosbach, K., and Schroder, U.; Preparation and Application of Magnetic Polymers for Targeting of Drugs, FEBS Letters, 112, (1979).
10. Newbower, R.; Magnetic Fluids in the Blood, IEEE Transactions on Magnetics, MAG.-9, 447, (1973).
11. Ranney, D.F., Senyei, A.E. and Widder, K.; Magnetically Responsive Microspheres and Other Carriers for the Biophysical Targeting of Antitumour Agents, Advances in Pharmacological Chemotherapy, 16; 213 (1979) (190 REF.).
12. Rembaum, A., and Dreyer, W.J.; Immunomicrospheres: Reagents for Cell Labeling and Separation, Science, 208, 364 (1980).
13. Scarpelli, D., Senyei, A. and Widder, K.; Magnetic Microspheres: A Model System for Site Specific Drug Delivery in Vivo, Proceedings of the Society for Experimental Biology and Medicine, 58, 141 (1978).
14. Senyei, A.E. and Widder, K.; Drug Targeting: Magnetically Responsive Albumin Microspheres-A Review of the System to Date; Gynecology and Oncology, 12 (1): (1981).
15. Senyei, AE., and Widder, K.; Magnetic Microspheres, 1. of Histochemical Cytochemistry, (1981).
16. Sugibayashi, K.; Biomedical Applications of Magnetic Fluids, Biomaterials, (1982).
17. Margel, S., Beitler, U., and Ofarim, M.; Polyacrolein Microspheres As A New Tool In Cell Biology, J. Cell Sci. 56, 157 (1982).

For Investigational Use Only. The performance characteristics of this product have not been established. This product is not sold for use in products for which prolonged contact with skin or implantation in the human body. Ferrotec Corporation and its group companies do not recommend this material as safe and effective for such uses and assumes no liability for any such use. Ferrotec Corporation and its group companies also do not warrant and

specifically disclaims that the use of its ferrofluid & magnetic particles for biomedical applications covered by another party's patent do not or will not infringe the intellectual property rights of the owner of such biomedical patent.

For further information about this product, contact Ferrotec or distributors in your area.

Ferrotec Corporation (JAPAN)

Tel: +81 479 73-6752 Fax: +81 479 73-6602 Email: ffsales@ferrotec.co.jp

Ferrotec (USA) Corporation

Tel: +1 603 883-9800 Fax: +1 603 883-2308 Email: sales@ferrotec.com

Ferrotec GmbH

Tel: +49 7022 9270-201 Fax: +49 7022 9270-1299 Email: fluid@de.ferrotec.com